


Oriel® Sol3A™ Class AAA Solar Simulators



- Output beam sizes 2x2", 4x4", 6x6", 8x8", and 12 x 12"
- Factory certified Class AAA CW systems
- Calibration certificate validating Class AAA performance for all 3 standards: IEC, ASTM and JIS
- Long-lived, highly reliable instruments designed specifically for 24/7 production environments
- Integrated variable attenuator provides capability to vary output from 0.1 to 1.0 suns
- Easy lamp replacement
- Non-reflective black finish reduces stray light
- Temperature sensors and interlocks ensure operator safety
- Improved optical design for maximum spatial uniformity
- Improved working distances accommodate larger samples
- Universal AC mains operation
- CE certified

Leveraging over 40 years of experience in light source and power supply design, Newport's Oriel® brand is proud to introduce the latest innovation in solar simulators – the Oriel Sol3A™ Class AAA Solar Simulator family. The product family includes:

Model 94123A = 12"x12" Beam Size

Model 94083A = 8"x8" Beam Size

Model 94063A = 6"x6" Beam Size ←

Model 94043A = 4"x4" Beam Size

Model 94023A = 2"x2" Beam Size

All Oriel Sol3A simulators are certified to IEC 60904-9 Edition 2 (2007), JIS C 8912, and ASTM E 927-05 standards for Spectral Match, Non-Uniformity of Irradiance, and Temporal Instability of Irradiance. By convention, Class AAA is reported with the first letter representing Spectral performance, the second letter Uniformity of Irradiance, and the third letter Temporal Stability. The Oriel Sol3A simulators all use a single lamp design to meet not one or two, but all three performance criteria without compromising the 1 Sun output power, providing true Class AAA performance. The Oriel Sol3A uses a black non-reflective finish to minimize stray light and incorporates captive screws for all panels requiring user access to facilitate lamp replacement, alignment, and filter changes. Safety interlocks prevent inadvertent exposure to UV light. The Oriel Sol3A rugged design is backed by the Newport Corporation's world wide organization.

WHY CLASS AAA?

Photovoltaic (PV) cell research and manufacturing are evolving at an amazing rate, with new technologies like thin film processes, tandem junction and multi-junction cells, organic thin films and dye sensitized cells all requiring careful evaluation of their performance. The old designation of "Class A" solar simulators had become standard in most test laboratories and manufacturing environments, but the definition of Class A has degraded to mean a system that is capable of meeting any one (typically spectral match) of the three performance parameters called out in the test methods. In October 2007 the IEC revised the method used to evaluate the performance of a solar simulator. The new standard calls for quantifying and reporting the performance of each of three performance criteria separately. A solar simulator is now measured as class A, B, or C for Spectral Match, Non-Uniformity of Irradiance, and Temporal Instability of Irradiance. The new Sol3A solar simulators have been designed to meet Class A performance for all 3 of the test requirements, Spectral Match, Non-Uniformity of Irradiance, and Temporal Instability of Irradiance as called out by IEC 60904-9 Edition 2 (2007) edition, JIS C 8912, and ASTM E 927-05.

By ensuring measurement uniformity that allows results comparability and traceability, Class AAA systems reduce binning variability of photovoltaic cell testing as compared to so called Class A, Class B or non-classified sources. This performance consistency allows for precise comparison of performance data for researchers engaged in developing novel solar materials.

Newport Corporation puts each Oriel Sol3A™ through rigorous testing for all 3 standards (IEC, STM, and JIS) to insure compliance and supplies a certificate of calibration for all three standards.

Oriel
INSTRUMENTS
A Newport Company

Oriel® Sol3A™ Class AAA Solar Simulators

DEFINING CLASS AAA PERFORMANCE STANDARDS

Photovoltaic standards mandate that Class AAA solar simulators meet demanding requirements in three key performance areas: spectral match to the solar spectrum, spatial non-uniformity of irradiance, and temporal instability of irradiance.

There are three standards that define solar simulator performance.

- IEC 60904-9 Edition 2 (2007) Photovoltaic Devices – Part 9: Solar Simulator Performance Requirements
- JIS C 8912-1998, Solar Simulators for Crystalline Solar Cells and Modules
- ASTM E 927-05 (2005) Specification for Solar Simulation for Terrestrial PV Testing

Table 1 Class AAA Standards and Specifications

Organization	IEC	JIS	ASTM
Performance Parameter	60904-9-2007	C 8912	E92-05
Spectral Match (fraction of ideal percentage)	0.75 - 1.25	0.75 - 1.25	0.75 - 1.25
Non-Uniformity of Irradiance	2.0%	≤±2%	2%
Temporal Instability	0.5% STI <2.0% LTI	≤±1%	2%

SPECTRAL MATCH

The standards define the spectral match of a solar simulator as a percentage of the integrated intensity in 6 spectral ranges (listed in Table 2). Any deviation from the specified percentages must then lie within a range that determines the class of the simulator. For Class AAA, this range is 0.75 to 1.25 times the ideal percentage.




Fig. 1 Oriel Sol3A Spectral Match with AM 1.5G spectral correction filter meets IEC, JIS, ASTM Class A requirements to for spectral match.

Spectral Match

Table 2 Ideal Spectral Match Defined by IEC Standards

Spectral Range (nm)	Total Irradiance Range (%)	Ideal %
400 - 500	13.9 - 23.1	18.5
500 - 600	15.1 - 25.1	20.1
600 - 700	13.7 - 22.9	18.3
700 - 800	11.1 - 18.5	14.8
800 - 900	9.2 - 15.3	12.2
900 - 1100	12.1 - 20.1	16.1




Fig. 2 Measured Uniformity of a 2x2" Oriel Sol3A Solar Simulator.

SPATIAL UNIFORMITY OF IRRADIANCE

The irradiance uniformity over the work area is the most difficult Class AAA requirement to achieve and maintain. Hot spots can lead to significant errors in measured cell efficiency and can cause inaccurate binning of cells. The Class AAA spatial non-uniformity performance standard is designed to minimize the impact of hot spots and has a very stringent requirement of ≤2%. The plot below shows the uniformity of the irradiance across a typical simulator working area. Each unit will come with a plot of irradiance non-uniformity. The working distance ranges for each simulator are listed in Table 3.

SOL3A WORKING DISTANCES

Table 3 Working Distances

Model	Size	Working Distance Range (Inches)
94023A	2 x 2 (50.8 x 50.8)	12.0 ±0.5
94043A	4 x 4 (101.6 x 101.6)	6.0 ±0.5
94063A	6 x 6 (152.4 x 152.4)	7.0 ±0.5
94083A	8 x 8 (203.2 x 203.2)	15.0 ±0.5
94123A	12x12 (304.8 x 304.8)	12.0 ±0.5
94123A-CPV	12x12 (304.8 x 304.8)	4.0 ±0.5




Fig. 3 Typical Output Variation of a 1.6 kW Oriel Sol3A Solar Simulator Over Time.

TEMPORAL INSTABILITY

Temporal Instability is the third performance parameter of Class AAA standards. It requires that the output light be stable over time in order to ensure that the lamp fluctuations do not distort the measurement of solar cell efficiency. Oriel's Sol3A easily meets the requirements for temporal instability as defined by the IEC, ASTM, and JIS standards. The IEC 60904-9 (2007) has the most stringent requirements for short term instability with a maximum allowable level of 0.5%. Figure 4 shows a typical instability response for an Oriel Sol3A utilizing a 50 millisecond data acquisition time.

Elapsed Time (sec)	60
No. of Data Points	1200
% Instability	0.369

ORIEL CLASS AAA SOLAR SIMULATOR KEY COMPONENTS

Illuminator Housing

The illuminator housing provides a safe enclosure for the lamp. Its powder coated flat black finish provides a durable surface that also minimizes stray light. It is equipped with safety interlock systems to ensure operator and system safety. Panels containing user serviceable components use captive screws that require no tools for ease of removal and replacement. Integral fan(s) provide forced air-cooling to maintain optimal lamp, optics and housing temperature. The Sol3A housing utilizes a design which allows the head orientation to produce a downward, sideways or upwards facing beam on the 2x2 or 4x4 models. Contact your Sales Engineer for details.

Integrated Shutter

The Oriel Sol3A Solar Simulator includes an upgraded shutter for production-environment operation. The newly designed shutter for the Class AAA systems is a rugged, single-blade shutter designed for >1 million cycles. Historically, our real-world performance has exceeded 10 million cycles on units in the field for many years. The shutter has a minimum exposure time of 200 ms and can be controlled via a contact closure or logic level input, or a convenient push-button switch on the illuminator housing.

Xenon Arc Lamp

The Oriel Sol3A Solar Simulator source is a CW system. This enables testing of all cell materials unlike flash-lamp based systems that are limited by the response time of the material allowing the cell to be soaked at a constant light level prior to testing. The lamp is an ozone-free xenon short arc lamp. We certify each source with the supplied lamp. For continuous production environments, we suggest purchasing replacement lamps and an alignment palette to ensure Class AAA compliance as lamps are replaced.

Air Mass 1.5G Filter

The combination of lamp and air mass filter produces the characteristic Class AAA spectra. Our Air Mass 1.5G Filter retains its optical properties under the conditions encountered without degradation of the filter.

POWER SUPPLY SPECIFICATIONS

Model	69920	69920	69920	69922	69922
Power Factor	>0.99	>0.99	>0.99	>0.99	>0.99
Input Voltage	95 - 264 VAC	95 - 264 VAC	95 - 264 VAC	190 - 264 VAC	190 - 264 VAC
Input Frequency	47 - 63 Hz	47 - 63 Hz	47 - 63 Hz	47 - 63 Hz	47 - 63 Hz
Output Power	350 - 1200 W	350 - 1200 W	350 - 1200 W	800 - 1800 W	800 - 1800 W
Output Current	17.5 - 55 A	17.5 - 55 A	17.5 - 55 A	25 - 70 A	25 - 70 A
Output Voltage Range	0 - 45 VDC	0 - 45 VDC	0 - 45 VDC	0 - 30 VDC	0 - 30 VDC
Line Regulation	0.01 %	0.01 %	0.01 %	0.01 %	0.01 %
Output Voltage Ripple	< 0.1 % rms	< 0.1 % rms	< 0.1 % rms	< 0.1 % rms	< 0.1 % rms
Meter Accuracy (% of full scale)	< 0.05%	< 0.05%	< 0.05%	< 0.05%	< 0.05%
Digital Meter Resolution, Voltage	0.1 VDC	0.1 VDC	0.1 VDC	0.1 VDC	0.1 VDC
Digital Meter Resolution, Power	1 W	1 W	1 W	1 W	1 W
Digital Meter Resolution, Current	0.1 A	0.1 A	0.1 A	0.1 A	0.1 A
Safety Interlock Voltage	12 Vdc/GND	12 Vdc/GND	12 Vdc/GND	12 Vdc/GND	12 Vdc/GND
Operating Mode	Constant current or constant power	Constant current or constant power	Constant current or constant power	Constant current or constant power	Constant current or constant power
Ambient Operating Temperature	0 - 45 °C	0 - 45 °C	0 - 45 °C	0 - 45 °C	0 - 45 °C
Weight	20 (9) [lb (kg)]	20 (9) [lb (kg)]	20 (9) [lb (kg)]	20 (9) [lb (kg)]	20 (9) [lb (kg)]

1. Uniformity is defined as:

$$\text{Uniformity (\%)} = (\text{Max irradiance} - \text{Min irradiance}) / (\text{Max irradiance} + \text{Min irradiance}) \times 100\%$$

Uniformity is measured using following Methods:

- Compliant to IEC 60904-9 2007 edition - section 5.3, 64 points in an 8x8 grid with detector size no bigger than each grid size.
- Compliant to JIS C 8912 - section 4.2, 17 points in a pre-defined pattern with detector size no bigger than 2 x 2 cm square.
- Compliant to ASTM E 927-05 section 8.3, 64 points in an 8x8 grid with detector size no bigger than each grid size.

2. Temporal Instability is defined as:

$$\text{Temporal Instability (\%)} = (\text{Max irradiance} - \text{Min irradiance}) / (\text{Max irradiance} + \text{Min irradiance}) \times 100\%$$

Temporal instability is measured by taking 20 samples per second for 60 seconds.

3. Spectral match is measured by using a calibrated spectroradiometer with a grating monochromator and a discrete detector. Measurement resolution is set at 2 nm.

Ordering Information

SOL3A ORDERING INFORMATION


Model	Description
94023A	Sol3A Class AAA Solar Simulator IEC/JIS/ASTM, 450 Watt Xenon, 2x2
94043A	Sol3A Class AAA Solar Simulator IEC/JIS/ASTM, 450 Watt Xenon, 4x4
94063A	Sol3A Class AAA Solar Simulator IEC/JIS/ASTM, 1000 Watt Xenon, 6x6
94083A	Sol3A Class AAA Solar Simulator IEC/JIS/ASTM, 1600 Watt Xenon, 8x8
94123A	Sol3A Class AAA Solar Simulator IEC/JIS/ASTM, 1600 Watt Xenon, 12x12
94123A-CPV	Sol3A Class AAA Solar Simulator IEC/JIS/ASTM, 1600 Watt Xenon, 12x12, Concentrator
940X3A-P	Variable Aperture for partial sun illumination on Sol3A solar simulators

REPLACEMENT LAMPS AND FILTERS


Model	Description
6280NS	450 Watt Xenon Short Arc Lamp, Ozone Free (94023A or 94043A)
6272	1000 Watt Xenon Short Arc Lamp, Ozone Free (94063A)
62726	1600 Watt Xenon, Ozone Free Arc Lamp (94083A)
62726	1600 Watt Xenon, Ozone Free Arc Lamp (94123A)
81368	Air Mass Filter 1.5 Global (94123A)

Contact a Sales Engineer for pricing of optional accessories.


Oriel® Sol3A™ Class AAA Solar Simulators




Dimensional diagram of 94023A Sol3A Solar Simulator.



Dimensional diagram of 94043A Sol3A Solar Simulator.




Dimensional diagram of 94063A Sol3A Solar Simulator.



Dimensional diagram of 94083A Sol3A Solar Simulator.

Oriel® Sol3A™ Class AAA Solar Simulators



A Newport Company

Newport Corporation, Global Headquarters
1791 Deere Avenue, Irvine, CA 92606, USA

www.newport.com

PHONE: 1-800-222-6440 1-949-863-3144 FAX: 1-949-253-1680 EMAIL: sales@newport.com
Complete listings for all global office locations are available online at www.newport.com/contact

	PHONE	EMAIL		PHONE	EMAIL
Belgium	+32-(0)0800-11 257	belgium@newport.com	Irvine, CA, USA	+1-800-222-6440	sales@newport.com
China	+86-10-6267-0065	china@newport.com	Netherlands	+31-(0)30 6592111	netherlands@newport.com
France	+33-(0)1-60-91-68-68	france@newport.com	United Kingdom	+44-1235-432-710	uk@newport.com
Japan	+81-3-3794-5511	spectra-physics@splasers.co.jp	Germany / Austria / Switzerland	+49-(0)6151-708-0	germany@newport.com
Taiwan	+886 -(0)2-2508-4977	sales@newport.com.tw			

Newport Corporation, Irvine and Santa Clara, California and Franklin, Massachusetts; Evry and Beaune-La-Rolande, France; Stahnsdorf, Germany and Wuxi, China have all been certified compliant with ISO 9001 by the British Standards Institution.

© 2012 Newport Corporation. All rights reserved. Oriel and the Newport logo are registered trademarks of Newport Corporation. Sol3A is a trademark of Newport Corporation.

DS-12082 (10/16)